

APP

产品信息

产品名称	产品货号	产品规格
APP	ORS0034PAb	100 μ l/50 μ l

产品简介

Amyloid precursor protein (APP), with 32-amino acid protein (about 35 kDa), is part of a super-family of transmembrane and secreted proteins. It appears to have a number of roles, including regulation of haemostasis and mediation of neuroprotection. APP undergoes alternative splicing, resulting in several isoforms. Proteolytic cleavage of amyloid precursor protein leads to the formation of the 4 kDa amyloid beta/A4 protein, which is present in human platelets. APP is also involved in the formation of neurofibrillary tangles and plaques that characterize the senile plaques of Alzheimer patients.

产品性质

Alias: APP

Clonality: Monoclonal Antibody

Antibody Type: Primary antibody

Swissprot: P05067

Clone: 1B11F3

Entrez GeneID: 351

Host\Isotype: Mouse IgG2a

Reactivity: Human

Purification: affinity purified by Protein G

Applications: IHC=1/200 - 1/1000; ELISA=1/10000

Immunogen: Purified recombinant fragment of APP expressed in E. Coli.

Formulation: Ascitic fluid containing 0.03% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信