

CBX2

产品信息

产品名称	产品货号	产品规格
CBX2	ORS0994PAb	100 μ l/50 μ l

产品简介

This gene encodes a component of the polycomb multiprotein complex, which is required to maintain the transcriptionally repressive state of many genes throughout development via chromatin remodeling and modification of histones. Disruption of this gene in mice results in male-to-female gonadal sex reversal. Mutations in this gene are also associated with gonadal dysgenesis in humans. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

产品性质

Alias: M33; CDCA6; SRXY5

Clonality: Monoclonal Antibody

Antibody Type: Primary antibody

Swissprot: Q14781

Clone: 4C11B10

Entrez GeneID: 84733

Host\Isotype: Mouse IgG1

Reactivity: Human, Mouse, Rat, Rabbit, Monkey

Purification: affinity purified by Protein G

Applications: WB=1/500 - 1/2000; IHC=1/200 - 1/1000; ICC=1/200 - 1/1000; FCM=1/200 - 1/400; ELISA=1/10000

Immunogen: Purified recombinant fragment of human CBX2 (AA: 402-525) expressed in E. Coli.

Formulation: Purified antibody in PBS with 0.05% sodium azide

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
 021-57875552

发票/服务/监督投诉电话: 13122716991



2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信