

## CD87

### 产品信息

产品名称	产品货号	产品规格
CD87	ORS7257PAb	50 $\mu$ l/100 $\mu$ l/1ml

### 产品简介

CD87, known as urokinase plasminogen activator-receptor (uPAR), is a 36-68 kD single-chain glycoposphotidylinositol (GPI)-anchored glycoprotein. The extracellular region of CD87 consists of three homologous domains, D1, D2, and D3. Membrane-bound and soluble (D2D3) forms of uPAR have been reported. The membrane-bound CD87 is expressed on granulocytes, monocytes/macrophages, dendritic cells, fibroblasts, endothelial cells and keratinocytes. CD87 plays a role in cell adhesion and chemotaxis through binding to  $\beta$ 1,  $\beta$ 2 and  $\beta$ 3 integrins. CD87-bound uPA converts plasminogen to plasmin.

### 产品性质

**Alias:** Urokinase plasminogen activator-receptor, uPA-R, uPAR

**Clonality:** Recombinant Mouse Monoclonal Antibody

**Antibody Type:** Primary antibody

**Swissprot:** Q03405

**Clone:** 8E2

**Entrez GeneID:** 5329

**Host\Isotype:** Mouse IgG2a

**Reactivity:** Human

**Purification:** affinity purified by Protein G

**Applications:** FCM=1/100-1/200

**Immunogen:** -

**Formulation:** Purified antibody in PBS with 0.05% sodium azide.

**Traits:** liquid

**Concentration:** -

**Storage:** Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

### 注意事项:

1. This product is for scientific research purposes only.



网站: <https://www.aoruisai.com>  
<https://www.aoruicell.com>



手机号: 13061929181(微信)  
13061929770(微信)



邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)



电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信