

Flag Tag

产品信息

产品名称	产品货号	产品规格
Flag Tag	ORS7270PAb	100 μ l/50 μ l

产品简介

—

产品性质

Alias: Flag;DYKDDDDK;DYKDDDDK Tag;

Clonality: Recombinant Mouse Monoclonal Antibody

Antibody Type: Tag Antibody

Swissprot: -

Clone: -

Entrez GeneID: -

Host\Isotype: -

Reactivity: Tag

Purification: affinity purified by Protein G

Applications: WB=1/5000-1/10000;IP=1/100-1/500;ICC=1/500-1/1000

Immunogen: Synthetic peptide DYKDDDDK (KLH-coupled)

Formulation: Purified antibody in TBS with 0.03% Proclin 300, 1%BSA and 50% glycerol.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aorucell.com> 13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信