

GST Tag

产品信息

产品名称	产品货号	产品规格
GST Tag	ORS0019PAb	100 μ l/50 μ l

产品简介

Epitope tagging offers an easy and universal strategy for the identification and purification of proteins derived by recombinant DNA technology. Expression vectors with Glutathione S-transferase (GST) tag have been developed widely for both prokaryotic and eukaryotic systems, because GST can provide both an easily detectable Tag and a simple purification process with little effect on the biological function of the protein encoded by inserted cDNA sequences. Antibodies of anti-GST tag are useful for checking protein expression using immunoblotting, immunoprecipitation and immunostaining techniques, as well as for immunoaffinity purification of protein expressed in vitro.

产品性质

Alias: -

Clonality: Monoclonal Antibody

Antibody Type: Tag Antibody

Swissprot: -

Clone: 6G9C6

Entrez GeneID: -

Host\Isotype: Mouse IgG1

Reactivity: Tag

Purification: affinity purified by Protein G

Applications: WB=1/500 - 1/2000; IP=1/10-1/20; IHC=1/100 - 1/1000; ICC=1/100 - 1/1000; FCM=1/200 - 1/400; ELISA=1/10000

Immunogen: Purified fusion protein with GST protein tag expressed in E. Coli.

Formulation: Purified antibody in PBS with 0.05% sodium azide

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:



网站: <https://www.aoruisai.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aorucell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信