

## His Tag

### 产品信息

产品名称	产品货号	产品规格
His Tag	ORS2402PAb	100 $\mu$ l/50 $\mu$ l

### 产品简介

The 6 X histidine motif is used as a tag on many recombinant proteins to facilitate purification. It can then be purified using immobilized metal ions, traditionally nickel, in a procedure termed immobilized-metal affinity chromatography. Further analysis of the protein can be achieved using an antibody specific to the 6 X histidine in techniques such as western blotting. The antibody recognizes the His-tag fused to targeted proteins in transfected or transformed cells.

### 产品性质

**Alias:** 6 His tag; HHHHHH tag

**Clonality:** Monoclonal Antibody

**Antibody Type:** Tag Antibody

**Swissprot:** -

**Clone:** 6A8H8

**Entrez GeneID:** -

**Host\Isotype:** Mouse IgG1

**Reactivity:** Tag

**Purification:** affinity purified by Protein G

**Applications:** WB=1/500-1/1000; IP=1/20-1/200; IHC=1/100 - 1/1000; ICC=1/100 - 1/1000; FCM=1/200 - 1/400; ELISA=1/10000

**Immunogen:** Synthesized peptide of His Tag (AA: HHHHHHHH-KLH).

**Formulation:** Purified antibody in PBS with 0.05% sodium azide

**Traits:** liquid

**Concentration:** -

**Storage:** Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

### 注意事项:

1. This product is for scientific research purposes only.



网站: <https://www.aoruisai.com>  
<https://www.aoruicell.com>



手机号: 13061929181(微信)  
 13061929770(微信)



邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)



电话: 021-57896695  
 021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信