

His Tag

产品信息

产品名称	产品货号	产品规格
His Tag	ORS2403PAb	100 μ l/50 μ l

产品简介

The 6 X histidine motif is used as a tag on many recombinant proteins to facilitate purification. It can then be purified using immobilized metal ions, traditionally nickel, in a procedure termed immobilized-metal affinity chromatography. Further analysis of the protein can be achieved using an antibody specific to the 6 X histidine in techniques such as western blotting. The antibody recognizes the His-tag fused to targeted proteins in transfected or transformed cells.

产品性质

Alias: 6 His tag; HHHHHH tag

Clonality: Monoclonal Antibody

Antibody Type: Tag Antibody

Swissprot: -

Clone: 2G7D9

Entrez GeneID: -

Host\Isotype: Mouse IgG1

Reactivity: Tag

Purification: affinity purified by Protein G

Applications: WB=1/500-1/1000; IP=1/10-1/20; IHC=1/100 - 1/1000; ICC=1/100 - 1/1000; FCM=1/200 - 1/400; ELISA=1/10000

Immunogen: Synthesized peptide of His Tag (AA: HHHHHHHH-KLH).

Formulation: Purified antibody in PBS with 0.05% sodium azide

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.



网站: <https://www.aoruisai.com> <https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信