

## Human Ig kappa chain - APC

### 产品信息

产品名称	产品货号	产品规格
Human Ig kappa chain - APC	ORS7125PAb	100 $\mu$ l/500 $\mu$ l/1ml

### 产品简介

Mouse Monoclonal Secondary Antibody to Human Ig kappa chain APC labeled, The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads., By immunoelectrophoresis and ELISA this antibody reacts specifically with Human Ig kappa chain. No antibody was detected against non immunoglobulin serum proteins.

### 产品性质

**Alias:** Mouse Anti-Human Ig kappa chain - APC; Mouse Anti Human Ig kappa chain APC

**Clonality:** Monoclonal Antibody

**Antibody Type:** Secondary Antibody

**Swissprot:** -

**Clone:** -

**Entrez GeneID:** -

**Host\Isotype:** Mouse

**Reactivity:** Human

**Purification:** affinity purified by Protein G

**Applications:** IHC=1/100 - 1/1000; ICC=1/100 - 1/1000; FCM=1/1000 - 1/4000; ELISA=Use at an assay dependent concentration

**Immunogen:** Human Ig kappa chain

**Formulation:** Purified antibody in PBS with 0.05% sodium azide, 1% BSA and 50% glycerol.

**Traits:** liquid

**Concentration:** -

**Storage:** Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

### 注意事项:

1. This product is for scientific research purposes only.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信