

Human Ig lambda chain - Cy5.5

产品信息

产品名称	产品货号	产品规格
Human Ig lambda chain - Cy5.5	ORS7108PAb	100 μ l/500 μ l/1ml

产品简介

Mouse Monoclonal Secondary Antibody to Human Ig lambda chain Cy5.5 labeled, The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads., By immunoelectrophoresis and ELISA this antibody reacts specifically with Human Ig lambda chain. No antibody was detected against non immunoglobulin serum proteins.

产品性质

Alias: Mouse Anti-Human Ig lambda chain - Cy5.5; Mouse Anti Human Ig lambda chain Cy5.5

Clonality: Monoclonal Antibody

Antibody Type: Secondary Antibody

Swissprot: -

Clone: -

Entrez GeneID: -

Host\Isotype: Mouse

Reactivity: Human

Purification: affinity purified by Protein G

Applications: IHC=1/100 - 1/1000; ICC=1/100 - 1/1000; FCM=1/1000 - 1/4000; ELISA=Use at an assay dependent concentration

Immunogen: Human Ig lambda chain

Formulation: Purified antibody in PBS with 0.05% sodium azide, 1% BSA and 50% glycerol.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信