

Human IgG Fc - Cy7

产品信息

产品名称	产品货号	产品规格
Human IgG Fc - Cy7	ORS7111PAb	100 μ l/500 μ l/1ml

产品简介

Mouse Monoclonal Secondary Antibody to Human IgG Fc Cy7 labeled, The antibody was isolated from ascitic by immunoaffinity chromatography using antigens coupled to agarose beads., By immunoelectrophoresis and ELISA this antibody reacts specifically with Human IgG Fc. No antibody was detected against non immunoglobulin serum proteins.

产品性质

Alias: Mouse Anti-Human IgG Fc - Cy7; Mouse Anti Human IgG Fc Cy7

Clonality: Monoclonal Antibody

Antibody Type: Secondary Antibody

Swissprot: -

Clone: -

Entrez GeneID: -

Host\Isotype: Mouse

Reactivity: Human

Purification: affinity purified by Protein G

Applications: IHC=1/100 - 1/1000; ICC=1/100 - 1/1000; FCM=1/1000 - 1/4000; ELISA=Use at an assay dependent concentration

Immunogen: Human IgG Fc

Formulation: Purified antibody in PBS with 0.05% sodium azide, 1% BSA and 50% glycerol.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信