

IL13

产品信息

产品名称	产品货号	产品规格
IL13	ORS7262PAb	50 μ l/100 μ l/1ml

产品简介

IL-13 is an immunoregulatory cytokine produced primarily by activated Th2 lymphocytes. IL-13 shares 30% amino acid sequence homology with IL-4 and demonstrates similar biological activities. The biological activities of IL-13 include: suppression of macrophage cytotoxic activity, upregulation of IL-1RA expression, and suppression of proinflammatory cytokine secretion.

产品性质

Alias: Interleukin-13, NC30

Clonality: Recombinant Mouse Monoclonal Antibody

Antibody Type: Primary antibody

Swissprot: P35225

Clone: 5E4

Entrez GeneID: 3596

Host\Isotype: Mouse IgG2a

Reactivity: Human

Purification: affinity purified by Protein G

Applications: FCM=1/100-1/200

Immunogen: -

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信