

IL17a

产品信息

产品名称	产品货号	产品规格
IL17a	ORS7263PAb	50 μ l/100 μ l/1ml

产品简介

IL-17A was initially identified from a subtracted cDNA library between closely related murine lymphoid cells called CTLA-8. IL-17A shares 58% homology with an open reading frame of the T-lymphotropic Herpes virus Samirii virus (viral IL-17). IL-17A belongs to a family of cytokines, which has five members: designated IL-17A-F. IL-17 is expressed by a unique lineage of CD4 T cells known as Th17 cells. This lineage selectively produces proinflammatory cytokines including IL-17, IL-17F, IL-21, and IL-22. In mice, the differentiation of this new lineage is initiated by TGF β -1 and IL-6 or IL-21, which activates Stat3 and induces the expression of the transcription factor retinoic acid-related orphan receptors (ROR γ t and ROR α). IL-23, which also activates Stat3, serves to maintain Th17 cells lineage.

产品性质

Alias: CTLA-8, CTLA8, IL-17, IL-17A, IL17, ILA17

Clonality: Recombinant Mouse Monoclonal Antibody

Antibody Type: Primary antibody

Swissprot: Q16552

Clone: 3F8

Entrez GeneID: 3605

Host\Isotype: Mouse IgG2a

Reactivity: Human

Purification: affinity purified by Protein G

Applications: FCM=1/100-1/200

Immunogen: -

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.



网站: <https://www.aoruisai.com> <https://www.aorucell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信