

LMNB1

产品信息

产品名称	产品货号	产品规格
LMNB1	ORS2129PAb	100 μ l/50 μ l

产品简介

<p>This gene encodes one of the two B-type lamin proteins and is a component of the nuclear lamina. A duplication of this gene is associated with autosomal dominant adult-onset leukodystrophy (ADLD). Alternative splicing results in multiple transcript variants. </p>

产品性质

Alias: LMN; ADLD; LMN2; LMNB; MCPH26

Clonality: Monoclonal Antibody

Antibody Type: Loading Control Antibody

Swissprot: P20700

Clone: 2A11A2

Entrez GeneID: 4001

Host\Isotype: Mouse IgG1

Reactivity: Human, Mouse, Rat

Purification: affinity purified by Protein G

Applications: WB=1/500 - 1/2000; IHC=1/200 - 1/1000; ICC=1/200 - 1/1000; FCM=1/200 - 1/400; ELISA=1/10000

Immunogen: Purified recombinant fragment of human LMNB1 (AA: 413-583) expressed in HEK293-6e cells supernatant.

Formulation: Purified antibody in PBS with 0.05% sodium azide

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信