

MMP14

产品信息

产品名称	产品货号	产品规格
MMP14	ORS0983PAb	100 μ l/50 μ l

产品简介

Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP's are secreted as inactive proproteins which are activated when cleaved by extracellular proteinases. However, the protein encoded by this gene is a member of the membrane-type MMP (MT-MMP) subfamily; each member of this subfamily contains a potential transmembrane domain suggesting that these proteins are expressed at the cell surface rather than secreted. This protein activates MMP2 protein, and this activity may be involved in tumor invasion.

产品性质

Alias: MMP-14; MMP-X1; MT-MMP; MT1MMP; MTMMP1; WNCHRS; MT1-MMP; MT-MMP 1

Clonality: Monoclonal Antibody

Antibody Type: Primary antibody

Swissprot: P50281

Clone: 6A11D1

Entrez GeneID: 4323

Host\Isotype: Mouse IgG1

Reactivity: Human

Purification: affinity purified by Protein G

Applications: IHC=1/200 - 1/1000;FCM=1/200 - 1/400;ELISA=1/10000

Immunogen: Purified recombinant fragment of human MMP14 (AA: 112-246) expressed in E. Coli.

Formulation: Purified antibody in PBS with 0.05% sodium azide

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.



网站: <https://www.aoruisai.com>
<https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信