

MSLN

产品信息

产品名称	产品货号	产品规格
MSLN	ORS7267PAb	50 μ l/100 μ l/1ml

产品简介

Human mesothelin, also known as CAK1 and MSLN, is a cell-surface protein that is expressed on normal mesothelial cells lining the pleura, pericardium, and peritoneum. It is expressed as a 70 kD pre-cursor protein with a furin cleavage site, which produces two products: 1) megakaryocyte potentiating factor (MPF) (30 kD) and 2) mesothelin (40 kD). MPF is a cytokine that is released. Mesothelin, for the most part, is membrane-bound by a glycosylphosphatidylinositol linkage and has three glycosylation sites. Soluble mesothelin is believed to result from proteolytic cleavage to release it from the membrane linkage, or from alternative splicing of the mesothelin gene.

产品性质

Alias: MPF, SMRP

Clonality: Recombinant Mouse Monoclonal Antibody

Antibody Type: Primary antibody

Swissprot: Q13421

Clone: 8G3

Entrez GeneID: 10232

Host\Isotype: Mouse IgG2a

Reactivity: Human

Purification: affinity purified by Protein G

Applications: FCM=1/100-1/200

Immunogen: -

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.



网站: <https://www.aoruisai.com> <https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信