

## TUBB3

### 产品信息

产品名称	产品货号	产品规格
<b>TUBB3</b>	<b>ORS0212PAb</b>	<b>100 <math>\mu</math> l/50 <math>\mu</math> l</b>

### 产品简介

Tubulin, beta 3, also known as TUBB3. Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non exchangeable site on the alpha-chain. Tubulin is a highly conserved protein with a molecular weight of ~50 kD. Microtubules play key roles in chromosome segregation in mitosis, intracellular transport, ciliary and flagellar bending, and structural support of the cytoskeleton. The two main tubulin isoforms,  $\alpha$ - and

### 产品性质

**Alias:** tubulin, beta 3; MC1R; TUBB4

**Clonality:** Monoclonal Antibody

**Antibody Type:** Loading Control Antibody

**Swissprot:** Q13509

**Clone:** 2E9

**Entrez GeneID:** 10381

**Host\Isotype:** Mouse IgG1

**Reactivity:** Human, Mouse, Rat

**Purification:** affinity purified by Protein G

**Applications:** WB=1/500-1/2000; ICC=1/200 - 1/1000; FCM=1/200 - 1/400; ELISA=1/10000

**Immunogen:** Purified recombinant fragment of human TUBB3 expressed in E. Coli.

**Formulation:** Ascitic fluid containing 0.03% sodium azide.

**Traits:** liquid

**Concentration:** -

**Storage:** Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

### 注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
 021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信