

Thioredoxin (TRX)

产品信息

产品名称	产品货号	产品规格
Thioredoxin (TRX)	ORS0075PAb	100 μ l/50 μ l

产品简介

Thioredoxin (TRX) is a small ubiquitous protein (MW12kDa) which is exist in a wide variety of prokaryotic and eukaryotic cells. Trx contains a redox active disulfide/dithiol group within the conserved Cys-Gly-Pro-Cys active site. This antibody is suitable for detecting fusion proteins which encode a Trx-Tag by immunoblotting and immunoprecipitation. The Monoclonal Antibody can detect a little Trx-Tag fusion proteins with negligible cross-reactivity with bacterial, insect, or mammalian lysates.

产品性质

Alias: Thioredoxin

Clonality: Monoclonal Antibody

Antibody Type: Tag Antibody

Swissprot: -

Clone: 1H6H6

Entrez GeneID: -

Host\Isotype: Mouse IgG1

Reactivity: Tag

Purification: affinity purified by Protein G

Applications: WB=1/500 - 1/2000; IHC=1/100 - 1/1000; ICC=1/100-1/200; FCM=1/200 - 1/400; ELISA=1/10000

Immunogen: Purified recombinant fusion protein with Thioredoxin (TRX) tag.

Formulation: Ascitic fluid containing 0.03% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.



网站: <https://www.aoruisai.com> <https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信