

ZAP70

产品信息

产品名称	产品货号	产品规格
ZAP70	ORS0033PAb	100 μ l/50 μ l

产品简介

ZAP70 (zeta-chain associated protein kinase), a 70 kDa member of the SYK tyrosine kinase family, plays a central role in lymphocyte activation and development, and is implicated in several immune disorders. ZAP70 controls TCR(T-cell antigen receptor)-linked signal transduction pathways. Its key role in thymocytes development and mature T lymphocytes activation has been illustrated by the characterization of several human immunodeficiencies presenting with mutations in the ZAP70 gene. ZAP70 is also expressed in several types of B-cell neoplasm and is easily detected by immunohistochemistry, providing a useful prognostic marker in patients with chronic lymphocytic leukaemia.

产品性质

Alias: SRK; STD; TZK

Clonality: Monoclonal Antibody

Antibody Type: Primary antibody

Swissprot: P43403

Clone: 3D2A4

Entrez GeneID: 7535

Host\Isotype: Mouse IgG1

Reactivity: Human

Purification: affinity purified by Protein G

Applications: WB=1/500 - 1/2000;IHC=1/200 - 1/1000;ELISA=1/10000

Immunogen: Purified recombinant fragment of human ZAP70 expressed in E. Coli.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.



网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信