

Olfactory Marker Protein

产品信息

产品名称	产品货号	产品规格
Olfactory Marker Protein	ORS5901PAb	100 μ l/50 μ l

产品简介

May act as a modulator of the olfactory signal-transduction cascade.

产品性质

Alias: Olfactory marker protein

Clonality: Monoclonal Antibody

Antibody Type: Primary antibody

Swissprot: Q64288

Clone: -

Entrez GeneID: 18378

Host\Isotype: Rabbit IgG

Reactivity: Mouse, Rat

Purification: affinity purified by Protein G

Applications: WB=1/500-1/1000; IP=1/10-1/20

Immunogen: Recombinant protein of mouse Olfactory Marker Protein

Formulation: Purified antibody in TBS with 0.05% sodium azide, 0.05% BSA and 50% glycerol.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

