

CAMSAP3 (N-term)

产品信息

产品名称	产品货号	产品规格
CAMSAP3 (N-term)	ORS10013MAb	100 μ l/50 μ l

产品简介

Microtubule minus-end binding protein that acts as a regulator of microtubule dynamics. Specifically required for zonula adherens biogenesis and maintenance by anchoring microtubules at their minus-ends to zonula adherens, leading to recruit KIFC3 kinesin to junctional site.

产品性质

Alias: Calmodulin-regulated spectrin-associated protein 3, Protein Nezha, CAMSAP3, KIAA1543

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: Q9P1Y5

Clone: -

Entrez GeneID: 57662

Host\Isotype: Rabbit IgG

Reactivity: Human, Mouse, Rat

Purification: affinity purified by Protein G

Applications: WB=1/1000

Immunogen: This CAMSAP3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 103-131 amino acids from the N-terminal region of human CAMSAP3.

Formulation: Purified antibody in PBS with 0.05% sodium azide.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信