

## METTL1

### 产品信息

产品名称	产品货号	产品规格
METTL1	ORS8733MAb	100 $\mu$ l/50 $\mu$ l

### 产品简介

This gene is similar in sequence to the *S. cerevisiae* YDL201w gene. The gene product contains a conserved S-adenosylmethionine-binding motif and is inactivated by phosphorylation.

### 产品性质

**Alias:** tRNA (guanine-N(7)-)-methyltransferase {ECO:0000255|HAMAP-Rule:MF\_03055}, 21133 {ECO:0000255|HAMAP-Rule:MF\_03055}, Methyltransferase-like protein 1 {ECO:0000255|HAMAP-Rule:MF\_03055}, tRNA (guanine(46)-N(7))-methyltransferase {ECO:0000255|HAMAP-Rule:MF\_030}

**Clonality:** Polyclonal Antibody

**Antibody Type:** Primary antibody

**Swissprot:** Q9UBP6

**Clone:** -

**Entrez GeneID:** 4234

**Host\Isotype:** Rabbit IgG

**Reactivity:** Human

**Purification:** affinity purified by Protein G

**Applications:** WB=1/1000

**Immunogen:** This METTL1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 147-175 amino acids from the Central region of human METTL1.

**Formulation:** Purified antibody in PBS with 0.05% sodium azide.

**Traits:** liquid

**Concentration:** -

**Storage:** Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

### 注意事项:

1. This product is for scientific research purposes only.



网站: <https://www.aoruisai.com>  手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)



邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)



电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)  
<https://www.aoruicell.com> 13061929770(微信)

邮箱: [sale4@aoruisai.com](mailto:sale4@aoruisai.com)  
[sale5@aoruisai.com](mailto:sale5@aoruisai.com)

电话: 021-57896695  
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信