

SERPINA3

产品信息

产品名称	产品货号	产品规格
SERPINA3	ORS0984MAb	100 μ l/50 μ l

产品简介

The protein encoded by this gene is a plasma protease inhibitor and member of the serine protease inhibitor class. Polymorphisms in this protein appear to be tissue specific and influence protease targeting. Variations in this protein's sequence have been implicated in Alzheimer's disease, and deficiency of this protein has been associated with liver disease. Mutations have been identified in patients with Parkinson disease and chronic obstructive pulmonary disease.

产品性质

Alias: ACT; AACT; GIG24; GIG25

Clonality: Polyclonal Antibody

Antibody Type: Primary antibody

Swissprot: P01011

Clone: -

Entrez GeneID: -

Host\Isotype: Rabbit IgG

Reactivity: Human

Purification: affinity purified by Protein G

Applications: WB=1/200-1/1000; IHC=1/25-1/50; ELISA=1/5000-1/10000

Immunogen: Synthetic peptide of human SERPINA3

Formulation: Purified antibody in PBS with 0.05% sodium azide and 50% glycerol.

Traits: liquid

Concentration: -

Storage: Store at 4° C short term. Aliquot and store at -20° C long term. Avoid freeze/thaw cycles.

注意事项:

1. This product is for scientific research purposes only.
2. Avoid repeated freeze-thaw cycles.
3. For your safety and health, please wear a lab coat and disposable gloves



网站: <https://www.aoruisai.com> <https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信