

重组蛋白

Catalog Number:ORSCF2434AMu

Recombinant *Bacillus subtilis* Uncharacterized protein

YuaB(yuaB) (Active)

Organism Species: *Bacillus subtilis* (strain 168)



使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: YUAB

Uniprot:P71014

Purity: Greater than 90% as determined by SDS-PAGE.

Endotoxin: Not test.

Source: in vitro E.coli expression system

Protein Length: Full Length

Expression Region: 1-181aa

Sequence:

MKRKLLSSLAISLSLGLLVSAPTASFAAESTSTKAHTESTMRTQSTASLFATITGASKT
EWSFSDIELTYRPNTLLSLGVMEFTLPSGFTANTKDTLNGNALRTTQILNNGKTVRVPLA
LDLLGAGEFKLKLNNKTLPAAGTYTFRAENKSLSIGNKFYAEASIDVAKRSTPPTQPCGC
N

Tag Info: N-terminal 6xHis-SUMO-tagged

Form:Lyophilized powder

Buffer formulation:Lyophilized from 10 mM Tris-HCl, 1 mM
EDTA, 6% Trehalose, pH 8.0

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:35.3 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

lyophilized form is 12 months at -20° C/-80° C.Avoid repeated

freeze/thaw cycles.



网站: <https://www.aoruisai.com>
<https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself



网站: <https://www.aoruisai.com>
<https://www.aoruicell.com>



手机号: 13061929181(微信)
13061929770(微信)



邮箱: sale4@aoruisai.com
sale5@aoruisai.com



电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信