

重组蛋白

Catalog Number:ORSEP0316R0

Recombinant Schmidtea mediterranea Activin 2, partial

Organism Species: Schmidtea mediterranea (Freshwater
planarian flatworm)

使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: N/A

Uniprot:U3RAB9

Purity: Greater than 85% as determined by SDS-PAGE.

Endotoxin: -

Source: E.coli

Protein Length: Partial

Expression Region: 64-406aa

Sequence:

NPVERKSSVYNDNAEEKIIQKDSFKDKVINDLKSKILQKLNKEAPKFNSSFDWSSMTEQ
IRSAALRLGVPFDETTENSQPTSESMIISLTKVNHCIDNSTNNQHGYDLTIPNMDNKILV
GANFIIQFLPLAFERTVEDNIIYNITINNLPEVYSFSVAQLKSNSLYIAYTNILKTLLD
NQRMSHSTFSFTVSISGDEQTTKLERMIDIENIVIELEYQVLPTPVRKRRDAGNTCIPNT
YFCCTQTLKVGIDELGWSQFIFAPTTFTINYCRGDCNSYLTFSN HARALNFARGRGLGS
NKDLRACCVPSYSSLTIMYLNENALEVSTFNDLIAATCGCM

Tag Info: N-terminal 6xHis-tagged

Form:Liquid or Lyophilized powder

Buffer formulation:-

Applications:-

Predicted Molecular Mass:-

Accurate Molecular Mass:42.9 kDa

Predicted isoelectric point:-

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20° C/-80° C, The shelf life of

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991




企业微信

企业微信

lyophilized form is 12 months at -20°C / -80°C . Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991

