

重组蛋白

Catalog Number:ORSEP5407AMu14

Recombinant Human Novel Coronavirus

Nucleoprotein(N) (D63G, R203M, D377Y) (Active)

Organism Species: Human Novel Coronavirus (SARS-CoV-2/
2019-nCoV)

使用说明手册V1.0

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信

【Product Information】

Gene Names: N

Uniprot: P0DTC9

Purity: Greater than 85.7% as determined by SDS-PAGE.

Endotoxin: Not test.

Source: E. coli

Protein Length: Full Length

Expression Region: 1-419aa (D63G, R203M, D377Y)

Sequence:

MSDNGPQNQRNAPRITFGGPSDSTGSNQNGERSGARSKQRRPQGLPNNTASWFTALTQHG
KEGLKFPRGQGVPINTNSSPDDQIGYYRRATTRIRGGDGKMKDLSRWYFYLLGTGPEAG
LPYGANKDGIWVATEGALNTPKDHIGTRNPANNAIIVLQLPQGTTLPKGFYAEGSRGGS
QASSRSSSRSRNSSRNSTPGSSMGTSPARMAGNGGDAALALLLDRLNQLESKMSGKGQQ
QQGQTVTKKSAAEASKKPRQKRTATKAYNVTQAFGRRGPEQTQGNFGDQELIRQGTDYKH
WPQIAQFAPSASAFFGMSRIGMEVTPSGTWLTYTGAIKLDDKDPNFKDQVILLNKHIDAY
KTFPPTEPKDKKKKKAYETQALPQRQKKQQTVTLLPAADLDDFSKQLQQSMSSADSTQA

Tag Info: N-terminal 6xHis-tagged

Form: Lyophilized powder

Buffer formulation: Lyophilized from 20 mM Tris-HCl, 0.5 M NaCl, 6% Trehalose, pH 8.0

Applications: -

Predicted Molecular Mass: -

Accurate Molecular Mass: 49.1 kDa

Predicted isoelectric point: -

【Storage and Stability】

Storage: Store at 2-8°C for one month, the shelf life of

liquid form is 6 months at -20°C/-80°C, The shelf life of

网站: <https://www.aoruisai.com> 手机号: 13061929181(微信)
<https://www.aoruicell.com> 13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991




企业微信

企业微信

lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. Avoid repeated freeze/thaw cycles.

【Identification】



【Important Note】

This product is for scientific research purposes only and should not be used for clinical diagnosis or other procedures, otherwise the consequences will be borne by oneself

网站: <https://www.aoruisai.com>
<https://www.aoruicell.com>

手机号: 13061929181(微信)
13061929770(微信)

邮箱: sale4@aoruisai.com
sale5@aoruisai.com

电话: 021-57896695
021-57875552

发票/服务/监督投诉电话: 13122716991



企业微信

企业微信